# CLEARFIL™ SILANE KIT

## ENGLISH INSTRUCTIONS FOR USE

### I. INTRODUCTION

CLEARFIL SILANE KIT is a surface treatment system for porcelain, hybrid ceramics and cured composite resin.

### II. INDICATIONS

CLEARFIL SILANE KIT is indicated for the following cases;

- [1] Surface treatment of prosthetic appliance made of porcelain, hybrid ceramics and cured composite resin.
- [2] Surface treatment of intraoral repair of fractured porcelain, hybrid ceramics or composite resin facing crowns/bridges using light-curing composite.

### III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers.

### IV. POSSIBLE SIDE EFFECTS

The mucous membrane may turn whitish when contacted by the CLEARFIL SE BOND Primer due to the coagulation of protein. This is a temporary phenomenon that will disappear in a few days.

### V. INCOMPATIBILITY

Do not use eugenol containing materials for pulp protection or temporary sealing, since the eugenol could retard the curing process.

### VI. PRECAUTIONS

#### 1. Safety precautions

- Avoid use of the product for patients with a history of hypersensitivity to methacrylate monomer.
- If any hypersensitivity occurs, such as a dermatitis, discontinue the use of the product and consult a physician.
- Use caution to prevent the product from coming in contact with the skin or getting into eye. Before using the product, cover the patient's eyes from splashing material.
- If the product comes in contact with human body tissue, take the following actions:
  - < If the product gets in the eye >  
Immediately wash the eye with copious amounts of water and consult a physician.
  - < If the product comes in contact with the skin >  
Immediately wipe it off with a cotton pledget moistened with alcohol or gauze and wash with copious amounts of water.
- Use caution to prevent the patient from accidentally swallowing the product.
- Avoid looking directly at the curing light when curing the product.
- Keep the product out of reach of children.

#### 2. Handling and manipulation precautions

- CLEARFIL SE BOND Primer polymerizes by visible light (especially ultraviolet rays). Use the light blocking plate to avoid an operating light or natural light (sunlight from windows) and use within 3 minutes from dispensing.
- The emitting tip of the visible light-curing activator should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
- Low intensity of light causes poor adhesion. Check the lamp for service life and light guide tip for contamination. It is advisable to check curing light intensity using an appropriate light-evaluating device at periodic intervals.

#### 3. Storage precautions

- Do not use after the expiration date. Note the expiration date on the outside of package.
- The product should be refrigerated when not in use; bring it to room temperature before using.
- Keep away from extreme heat or direct sunlight.
- Replace the bottle cap as soon as possible after resin has been dispensed from the bottle.

### VII. COMPONENTS

This Kit consists of one bottle each of CLEARFIL SE BOND Primer, CLEARFIL PORCELAIN BOND ACTIVATOR and K-ETCHANT GEL.

- 1) CLEARFIL SE BOND Primer  
1.5ml X 1 bottle
- 2) CLEARFIL PORCELAIN BOND ACTIVATOR  
1ml X 1 bottle
- 3) K-ETCHANT GEL  
1.5ml X 1 bottle

### VIII. CLINICAL PROCEDURES

#### A. Treatment of prosthetic appliances

Silanating treatment on prosthetic appliances made of esthetic materials (porcelain, hybrid ceramics or composite resin) can be accomplished using CLEARFIL SE BOND Primer in conjunction with CLEARFIL PORCELAIN BOND ACTIVATOR.

##### A-1. Acid etching of adherent surface

Apply phosphoric acid (K-ETCHANT GEL) to the adherent surface and leave it in place for 5 seconds before washing and drying.

##### A-2. Silane treatment

- 1) Dispense one drop each of CLEARFIL SE BOND Primer and CLEARFIL PORCELAIN BOND ACTIVATOR into a well of the mixing dish and mix them.

#### [NOTE]

CLEARFIL SE BOND Primer and CLEARFIL PORCELAIN BOND ACTIVATOR should be mixed immediately before application.

#### [CAUTION]

CLEARFIL SE BOND Primer polymerizes by visible light (especially ultraviolet rays). CLEARFIL SE BOND Primer or the mixture will set to a gel if left under the operating light or natural light. Use the light blocking plate to avoid operating light or natural light (sunlight from windows) and use within 3 minutes from dispensing.

- 2) Apply mixture to the facing material and metal surface with a sponge or a disposable brush tip.
  - 3) Leave it in place for 5 seconds. Use caution not to allow saliva or exudate to contact the treated surface.
  - 4) After that the 5 seconds, evaporate the volatile ingredients with a mild oil-free air stream. The treated metal surface must be sufficiently dried; if not, adhesion might be seriously impaired.
- [CAUTION]  
Do not wash the treated surface. If the treated surface is contaminated with saliva, wash it with water, dry, or clean with alcohol, and treated with the mixture again.

##### A-3. Cementing

Cement the prosthetic appliance using PANAVIA cement according to the manufacturer's instructions.

#### B. Surface treatment of intraoral repair of fractured facing crown made of porcelain (PFM), hybrid ceramics, or composite resin using light-curing composite.

##### B-1. Preparation of fractured surface

- 1) Facing material surface  
Using a diamond point, remove a layer of the fractured surface, and place a bevel at the marginal area.
- 2) Metal surface  
Roughen the metal surface with a diamond point.

##### B-2. Acid etching of facing material surface

Apply phosphoric acid (K-ETCHANT GEL) to the facing material surface and leave it in place for 5 seconds before washing and drying. If the adherent surface extends to uncut enamel, apply phosphoric acid (K-ETCHANT GEL) and leave it in place for 10 seconds before washing and drying.

##### B-3. Precious metal surface

When precious metal is used, apply ALLOY PRIMER according to the manufacturer's instructions.

##### B-4. Silane treatment

Use the same procedure described in A-2.

#### [CAUTION]

If the adherent surface extends to the tooth, leave mixture in place for 20 seconds.

##### B-5. Bonding

- 1) Dispense the necessary amount of CLEARFIL SE BOND Bond into a well of the mixing dish.

#### [CAUTION]

CLEARFIL SE BOND Bond polymerizes by visible light (especially ultraviolet rays). CLEARFIL SE BOND Bond will set to a gel if under the operating light or natural light. Use the light blocking plate to avoid operating light or natural light (sunlight from windows) and use within 3 minutes from dispensing.

- 2) Apply CLEARFIL SE BOND Bond to the entire fractured surface with a sponge or disposable brush tip.

- 3) Use a light air stream to make the bond film as uniform as possible.

- 4) Light-cure CLEARFIL SE BOND Bond for 10 seconds with a visible light-curing activator.

##### B-6. Light-cure opaquer applying

- 1) Use opaquer resin such as CLEARFIL ST OPAQUER for metal surface to prevent metal shine through.

- 2) Light-cure the composite for 40 seconds.

##### B-7. Light-cure composite filling

Apply composite (e.g. CLEARFIL AP-X) into the cavity, light-cure, finish and according to the manufacturer's instructions.

#### [NOTE]

Use opaquer resin such as CLEARFIL ST OPAQUER for metal surface to prevent metal shine through.

#### [WARRANTY]

KURARAY MEDICAL INC. will replace any product that is proved to be defective. KURARAY MEDICAL INC. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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CLEARFIL and PANAVIA are trademarks of KURARAY CO., LTD.
CLEARFIL y PANAVIA son marcas registradas de KURARAY CO., LTD.