CLEARFIL™ SE BOND

RESIN-BASED DENTAL ADHESIVE SYSTEM

ENGLISH INSTRUCTIONS FOR USE

I. INTRODUCTION

CLEARFIL SE BOND is a light-curing bonding system, and consists of a self-etching primer and a bonding agent. The primer offers simultaneous treatment of both dentine and enamel.

II. INDICATIONS

CLEARFIL SE BOND is indicated for the following cases:

- [1] Direct filling restorations using light-curing composite or compomer
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of hypersensitive and/or exposed root surfaces
- [4] Intraoral repairs of fractured facing crowns made of porcelain, hybrid ceramics or composite resin using light-curing composite
- [5] Surface treatment of prosthetic appliances made of porcelain, hybrid ceramics and cured composite resin

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers.

IV. POSSIBLE SIDE EFFECTS

The mucous membrane may turn whitish when contacted with PRIMER due to the coagulation of protein. This is a temporary phenomenon that will disappear in a few days.

V. INCOMPATIBILITY

Do not use eugenol containing materials for pulp protection and temporary sealing, as they retard the curing process.

VI. PRECAUTIONS

1. Safety precautions

1. Avoid use of the product for patients with a history of hypersensitivity to methacrylate monomers.
2. If any hypersensitivity such as dermatitis occurs, discontinue use of the product and consult a physician.
3. Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
4. If the product comes in contact with the human body, take the following measures:
<If the product gets in the eye>
Immediately wash the eye with copious amounts of water and consult a physician.
<If the product comes in contact with the skin>
Immediately wipe it off with a cotton pledget moistened with alcohol or gauze and wash with copious amounts of water.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. Avoid looking directly at the curing light when curing the product.
7. Keep the product out of reach of children.

2. Handling and manipulation precautions

1. PRIMER and BOND polymerize by visible light (especially ultraviolet rays). Use the light blocking plate to avoid operating light or natural light (sunlight from windows) and use within 3 minutes from dispensing.

2. The emitting tip of the visible light-curing activator should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
3. Low intensity of light causes poor adhesion. Check the lamp for service life and the visible light guide tip for contamination. It is advisable to check the light intensity of light-curing activator using an appropriate light-checker at periodic intervals.

3. Storage precautions

1. Do not use after the expiration date. Note expiration date on the outside of package.
2. The product should be refrigerated when not in use; bring it to room temperature before using.
3. Keep away from extreme heat or direct sunlight.
4. Replace the bottle cap as soon as possible after resin has been dispensed from the bottle.

VII. COMPONENTS

An introductory kit consists of a self-etching primer, a bonding agent and accessories.

1) PRIMER (self-etching primer) 6ml

Principle ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- 2-Hydroxyethyl methacrylate (HEMA)
- Hydrophilic dimethacrylate
- dl-Camphorquinone
- N,N-Diethanol-p-toluidine
- Water

2) BOND (bonding agent) 5ml

Principle ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- Bis-phenol A diglycidylmethacrylate (Bis-GMA)
- 2-Hydroxyethyl methacrylate (HEMA)
- Hydrophobic dimethacrylate
- dl-Camphorquinone
- N,N-Diethanol-p-toluidine
- Silanated colloidal silica

3) Accessories (in introductory kit)

- | | |
|-----------------------|---------------|
| Disposable brush tips | 50 pieces x 2 |
| Brush tip handles | 2 pieces |
| Mixing dish | 1 piece |
| Light blocking plate | 1 piece |
| Outer case | 1 piece |

* The PRIMER and BOND can be easily dispensed using this outer case. Only this case, which contains PRIMER and BOND, should be refrigerated when not in use.

VIII. CLINICAL PROCEDURES

- A. Direct filling restorations, cavity sealing and treatment of hypersensitive teeth or exposed root surface using light-curing composite (or compomer)
- A-1. Cleaning tooth structure
An insufficiently cleaned cavity may impair the adhesive performance. Make sure the cavity is sufficiently cleaned.
- A-2. Moisture control
In order to produce optimal results, avoid contamination of the treatment area from saliva or blood. A rubber dam is recommended to keep the tooth dry.
- A-3. Cavity preparations
Remove any infected dentine using CARIES DETECTOR and prepare the cavity in the usual manner.

A-4. Pulp protection

Any actual or near pulp exposure should be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

A-5. Acid etching uncut enamel

If there is a possibility of resin spreading over uncut enamel, apply K-ETCHANT GEL to the enamel, let it remain for 10 seconds, wash with water, and then dry.

[CAUTION]

The use of PRIMER alone does not sufficiently condition uncut enamel. Overfilling of resin onto unetched, uncut enamel could cause marginal discoloration.

A-6. Tooth surface treatment

1. Dispense the necessary amount of PRIMER into a well of the mixing dish immediately before application.

[CAUTION]

PRIMER polymerizes by visible light (especially ultraviolet rays). PRIMER will set to a gel if left under the natural light or operating light.

Use the light blocking plate to avoid operating light or natural light (sunlight from windows) and use within 3 minutes from dispensing.

2. Apply PRIMER to the entire cavity wall with a sponge or a disposable brush tip. Leave it in place for 20 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces for at least 20 seconds.
3. After conditioning the tooth surface for 20 seconds, evaporate the volatile ingredients with a mild oil-free air stream.

[CAUTION]

Avoid pooling of PRIMER. Also do not wash.

Observe the drying method, and treatment time to ensure optimum adhesion.

Avoid touching treated surface. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat with PRIMER again.

A-7. Bonding

1. Dispense the necessary amount of BOND into a well of the mixing dish.

[CAUTION]

BOND polymerizes by visible light (especially ultraviolet rays). BOND will set to a gel if left under the operating light or natural light.

Use the light blocking plate to avoid operating light or natural light (sunlight from windows) and use within 3 minutes from dispensing.

2. Apply BOND to the entire surface of the cavity with a sponge or a disposable brush tip.
3. After application, make the bond film as uniform as possible using a gentle oil-free air stream.

[CAUTION]

A strong air blast will scatter the bonding agent, resulting in poor adhesion.

4. Light-cure the BOND for 10 seconds with a visible light curing activator.

A-8a. Direct filling restorations using light-curing composite (or compomer)

Apply composite (e.g. CLEARFIL AP-X) or compomer into the cavity, light-cure, finish and polish according to the manufacturer's instructions for use.

A-8b. Cavity sealing and treatment of hypersensitive teeth or exposed root surfaces
Apply a thin coat of CLEARFIL LINER BOND 2 PROTECT LINER F to the tooth, and light-cure the composite for 20 seconds. Remove unpolymerized resin with a cotton pledget moistened with alcohol.

B. Intraoral repairs of fractured facing crowns made of porcelain (PFM), hybrid ceramics or composite resin using light-curing composite

B-1. Preparation of fractured surfaces

1) Facing material surface

Using a diamond point, remove a layer of the fractured surface, and place a bevel at the marginal area.

2) Metal surface

Roughen the metal surface with a diamond point.

B-2. Acid etching of facing material surface

Apply K-ETCHANT GEL to the facing material surface and leave it in place for 5 seconds before washing and drying. If the adherent surface extends to uncut enamel, apply K-ETCHANT GEL and leave it in place for 10 seconds before washing and drying.

B-3. Precious metal surface

When precious metal is used, apply ALLOY PRIMER according to the manufacturer's instructions for use.

B-4. Silane treatment

1. Dispense one drop each of PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR into a well of the mixing dish and mix them.

[NOTE]

PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR should be mixed immediately before application.

[CAUTION]

PRIMER polymerizes by visible light (especially ultraviolet rays). PRIMER or the mixture will set to a gel if left under the operating light or natural light.

Use the light blocking plate to avoid operating light or natural light (sunlight from windows) and use within 3 minutes from dispensing.

2. Apply mixture to the facing material and metal surfaces with a sponge or a disposable brush tip.

3. Leave it in place for 5 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.

4. If the adherent surface extends to the tooth, leave mixture in place for 20 seconds.

5. After applying PRIMER, evaporate the volatile ingredients with a mild oil-free air stream. The treated metal surface must be sufficiently dried; if not, adhesion might be seriously impaired.

[CAUTION]

Do not wash treated surface. If the treated surface is contaminated with saliva, wash it with water, dry, or clean with alcohol, and treat with the mixture again.

B-5. Bonding

1. Dispense the necessary amount of BOND into a well of the mixing dish.

[CAUTION]

BOND polymerizes by visible light (especially ultraviolet rays). BOND will set to a gel if left under the operating light or natural light.

Use the light blocking plate to avoid operating light or natural light (sunlight from windows) and use within 3 minutes from dispensing.

2. Apply BOND to the entire cavity wall with a sponge or a disposable brush tip.

3. Use a light air stream to make the bond film as uniform as possible.

4. Light-cure BOND for 10 seconds with a visible light-curing activator.

B-6. Light-curing composite filling

Apply composite (e.g. CLEARFIL AP-X) into the cavity, light-cure, finish and polish according to the manufacturer's instructions for use.

[NOTE]

Use opaque resin such as PHOTO CLEARFIL OPAQUER for metal surface to prevent metal shine through.

C. Treatment of prosthetic appliances

Silanating treatment on prosthetic appliances made of esthetic materials (porcelain, hybrid ceramics or composite resin) can be accomplished using PRIMER in conjunction with CLEARFIL PORCELAIN BOND ACTIVATOR.

C-1. Acid etching of adherent surface

Apply K-ETCHANT GEL to the adherent surface and leave it in place for 5 seconds before washing and drying.

C-2. Silane treatment

Use the same procedure described in B -4.

C-3. Cementing

Cement the prosthetic appliance using PANAVIA cement according to the manufacturer's instructions for use.

[CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

[WARRANTY]

KURARAY MEDICAL INC. will replace any product that is proved to be defective. KURARAY MEDICAL INC. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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CLEARFIL and PANAVIA are trademarks of KURARAY CO., LTD.
CLEARFIL y PANAVIA son marcas registradas de KURARAY CO., LTD.