# CLEARFIL™ PROTECT BOND

DENTAL BONDING AGENT

## ENGLISH INSTRUCTIONS FOR USE

### I. INTRODUCTION

CLEARFIL PROTECT BOND is a light-cure bonding system, consisting of a self-etching primer and a fluoride-releasing bonding agent. The primer simultaneously treats both dentin and enamel.

### II. INDICATIONS

CLEARFIL PROTECT BOND is indicated for the following uses:

- [1] Direct restorations using light-cured composite resin or compomer
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of hypersensitive and/or exposed root surfaces
- [4] Intraoral repairs of fractured crowns/bridges made of porcelain, hybrid ceramics or composite resin using light-cured composite resin
- [5] Surface treatment of prosthetic appliances made of porcelain, hybrid ceramics and cured composite resin
- [6] Core build-ups using light- or dual-cured composite resin
- [7] Cavity sealing under amalgam restorations

### III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers.

### IV. POSSIBLE SIDE EFFECTS

The mucous membrane may turn whitish when contacted by the PRIMER due to the coagulation of protein. This is a temporary phenomenon that will disappear in a few days.

### V. INCOMPATIBILITY

- [1] Do not use eugenol containing materials for pulp protection or temporary sealing, since the eugenol could retard the bonding system curing process.
- [2] Do not use hemostatics containing ferric compounds, since they may impair adhesion and the remaining ferric ion may cause discoloration in the layer or surrounding gingiva.

### VI. PRECAUTIONS

#### 1. Safety precautions

1. Avoid use of the product for patients with a history of hypersensitivity to methacrylate monomers.
2. If any hypersensitivity occurs, such as dermatitis, discontinue use of the product and consult a physician.
3. Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
4. If the product comes in contact with human body tissues, take the following actions:
  - <If the product gets in the eye>  
Immediately wash the eye with copious amounts of water and consult a physician.
  - <If the product comes in contact with the skin>  
Immediately wipe it off with a cotton pledget moistened with alcohol or gauze and wash with copious amounts of water.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. Avoid looking directly at the dental curing light when curing the product.
7. Keep the product out of reach of children.

#### 2. Handling and manipulation precautions

1. PRIMER and BOND will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes from dispensing.
2. The emitting tip of the dental curing light should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
3. Low intensity of light causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.

### 3. Storage precautions

1. Do not use after the expiration date. Note the expiration date on the outside of package.
2. The product shall be refrigerated when not in use; bring it to room temperature before using.
3. Keep away from extreme heat or direct sunlight.
4. Replace the bottle cap as soon as possible after resin has been dispensed from the bottle.

### VII. COMPONENTS

An introductory kit consists of a self-etching primer, a bonding agent and accessories.

- |                                 |     |
|---------------------------------|-----|
| 1) PRIMER (self-etching primer) | 6ml |
|---------------------------------|-----|

#### Principle ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- 12-Methacryloyloxydodecylpyridinium bromide (MDPB)
- 2-Hydroxyethyl methacrylate (HEMA)
- Hydrophilic dimethacrylate
- Water

- |                                  |     |
|----------------------------------|-----|
| 2) BOND (fluoride-bonding agent) | 5ml |
|----------------------------------|-----|

#### Principle ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- Bis-phenol A diglycidylmethacrylate (Bis-GMA)
- 2-Hydroxyethyl methacrylate (HEMA)
- Hydrophobic dimethacrylate
- dl-Camphorquinone
- N,N-Diethanol-p-toluidine
- Silanated colloidal silica
- Surface treated sodium fluoride

#### 3) Accessories (in introductory kit)

- |                         |               |
|-------------------------|---------------|
| · Disposable brush tips | 50 pieces x 2 |
| · Brush tip handles     | 2 pieces      |
| · Mixing dish           | 1 piece       |
| · Light blocking plate  | 1 piece       |
| · Outer case            | 1 piece       |

#### [NOTE]

- \* The PRIMER and BOND can easily be dispensed using the outer case container.
- \* The PRIMER and BOND shall be refrigerated when not in use.

### VIII. CLINICAL PROCEDURES

#### A. Direct restorations, cavity sealing and treatment of hypersensitive teeth or exposed root surfaces using light-cured composite resin (or compomer)

##### A-1. Cleaning tooth structure

An adequately cleaned cavity assures maximum adhesive performance. Be sure the cavity is adequately cleaned.

##### A-2. Moisture control

In order to produce optimal results, avoid contamination of the treatment area from saliva or blood. A rubber dam is recommended to keep the tooth clean and dry.

##### A-3. Cavity preparations

Remove any infected dentin using CARIES DETECTOR as a guide and prepare the cavity in the usual manner.

##### A-4. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

##### A-5. Acid etching uncut enamel

If there is a possibility of resin spreading over uncut enamel, apply K-ETCHANT GEL to the enamel, let it remain for 10 seconds, wash with water, and then dry.

#### [CAUTION]

The use of PRIMER alone may not sufficiently condition uncut enamel. Overfilling of resin onto unetched, uncut enamel could cause marginal discoloration.

##### A-6. Tooth surface treatment

1. Dispense the necessary amount of PRIMER into a well of the mixing dish immediately before application.

#### [CAUTION]

PRIMER will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes from dispensing.

2. Apply PRIMER to the entire cavity wall with a disposable brush tip. Leave it in place for 20 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces for at least 20 seconds.
3. After conditioning the tooth surface for 20 seconds, evaporate the volatile ingredients with a mild oil-free air stream.

#### [CAUTION]

Avoid pooling of PRIMER. Do not wash after applying PRIMER. Observe the drying method and treatment time to ensure optimum adhesion. Avoid touching the treated surface. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat with PRIMER again.

#### A-7. Bonding

1. Dispense the necessary amount of BOND into a well of the mixing dish.

#### [CAUTION]

BOND will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes from dispensing.

2. Apply BOND to the entire surface of the cavity with a disposable brush tip.
3. After applying BOND, create a uniform bond film using a gentle oil-free air flow.

#### [CAUTION]

A strong air blast will scatter the bonding agent, resulting in poor adhesion.

4. Light-cure the BOND for 10 seconds with a dental curing light (irradiate wave length:400-515nm).

#### A-8. Either A-8a or A-8b

##### A-8a. Direct restorations using light-cured composite resin (or compomer)

Apply composite resin (e.g. CLEARFIL AP-X) or compomer into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

##### A-8b. Cavity sealing and treatment of hypersensitive teeth or exposed root surfaces

Apply a thin coat of CLEARFIL LINER BOND 2 PROTECT LINER F to the tooth, and light-cure the composite resin for 20 seconds. Remove unpolymerized resin with a cotton pledget moistened with alcohol.

### B. Intraoral repairs of fractured crowns/bridges made of porcelain, hybrid ceramics or composite resin using light-cured composite resin

#### B-1. Preparation of fractured surfaces

- 1) Facing material surface  
Using a diamond point, remove a layer of the fractured surface, and place a bevel at the marginal area.
- 2) Metal surface  
Roughen the metal surface with a diamond point.

#### B-2. Acid etching of facing material surface

Apply K-ETCHANT GEL to the facing material surface and leave it in place for 5 seconds before washing and drying. If the adherent surface extends to uncut enamel, apply K-ETCHANT GEL to the enamel surface and leave it in place for 10 seconds before washing and drying.

#### B-3. Precious metal surface

When precious metal is used, apply ALLOY PRIMER according to the manufacturer's instructions.

#### B-4. Silane treatment

1. Dispense one drop each of PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR into a well of the mixing dish and mix them.

[NOTE]

PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR should be mixed immediately before application.

[CAUTION]

PRIMER or the mixture will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes from dispensing.

2. Apply mixture to the facing material and metal surfaces with a disposable brush tip.
3. Leave it in place for 5 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
4. If the adherent surface extends to the tooth, leave mixture in place for 20 seconds.
5. After applying PRIMER, evaporate the volatile ingredients with a mild oil-free air stream. The treated metal surface must be sufficiently dried; if not, adhesion might be seriously impaired.

[CAUTION]

Do not wash treated surface. If the treated surface is contaminated with saliva, wash it with water, dry, or clean with alcohol, and treat with the mixture again.

**B-5. Bonding**

1. Dispense the necessary amount of BOND into a well of the mixing dish.

[CAUTION]

BOND will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes from dispensing.

2. Apply BOND to the entire cavity wall with a disposable brush tip.
3. After applying BOND, create a uniform bond film using a gentle oil-free air flow.
4. Light-cure the BOND for 10 seconds with a dental curing light (irradiate wave length:400-515nm).

**B-6. Light-cured composite resin filling**

Apply composite resin (e.g. CLEARFIL AP-X) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

[NOTE]

Use opaque resin such as CLEARFIL ST OPAQUER for metal surface to prevent metal shine through.

**C. Treatment of prosthetic appliances**

Silanating treatment of prosthetic appliances made of esthetic materials (porcelain, hybrid ceramics or composite resin) can be accomplished using PRIMER in conjunction with CLEARFIL PORCELAIN BOND ACTIVATOR.

**C-1. Acid etching of adherent surface**

Apply K-ETCHANT GEL to the adherent surface and leave it in place for 5 seconds before washing and drying.

**C-2. Silane treatment**

Use the same procedure described in B-4.

**C-3. Cementing**

Cement the prosthetic appliance using PANAVIA cement according to the manufacturer's instructions.

**D. Core build-ups using light- or dual-cured composite resin**

**D-1. Cleaning tooth structure**

An adequately cleaned tooth surface assures maximum adhesive performance. Be sure the tooth surface is adequately cleaned.

**D-2. Moisture control**

In order to produce optimal results, avoid contamination of the treatment area from saliva or blood. A rubber dam is recommended to keep the tooth clean and dry.

**D-3. Tooth preparations**

Remove existing restorations and/or caries and place pin, post or matrix as needed.

**D-4. Pulp protection**

In case of pulp exposure or near pulp exposure, a hard setting calcium hydroxide material is recommended. Do not use eugenol materials for pulp protection.

**D-5. Application of primer**

1. Dispense the necessary amount of PRIMER into a well of the mixing dish immediately before application.

[CAUTION]

PRIMER will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes from dispensing.

2. Apply PRIMER to all the tooth surfaces (and the post if applicable) with a disposable brush tip. Leave it in place for 20 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces for at least 20 seconds.
3. After conditioning all the tooth surfaces (and the post) for 20 seconds, evaporate the volatile ingredients with a mild oil-free air stream.

[CAUTION]

Avoid pooling of PRIMER. Do not wash after applying PRIMER. Observe the drying method and treatment time to ensure optimum adhesion. Avoid touching the treated surface. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat with PRIMER again.

**D-6. Bonding**

1. Dispense the necessary amount of BOND into a well of the mixing dish.

[CAUTION]

BOND will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes from dispensing.

2. Apply BOND to all the tooth surfaces with a disposable brush tip.
3. After applying BOND, create a uniform bond film using a gentle oil-free air flow.

[CAUTION]

A strong air blast will scatter the bonding agent, resulting in poor adhesion.

4. Light-cure the BOND for 10 seconds with a dental curing light (irradiate wave length:400-515nm).

**D-7. Core build-ups**

1. Place light- or dual-cured composite resin for core-build up, e.g. CLEARFIL PHOTO CORE, using a hand instrument.
2. Light-cure the core material according to the manufacturer's instructions.
3. After curing, check the surface hardness with an explorer. Initial crown preparation may begin immediately.

**E. Cavity sealing under amalgam restorations**

**E-1. Cleaning tooth structure**

**E-2. Moisture control**

**E-3. Cavity preparations**

Prepare the cavity in the normal manner for an amalgam restoration.

**E-4. Pulp protection**

**E-5. Acid etching uncut enamel**

**E-6. Tooth surface treatment**

**E-7. Bonding**

Use the same procedures described in A-1 to A-7 for the above steps.

**E-8. Amalgam filling**

Mix amalgam powder and mercury, condense the mixture into the cavity, finish and polish according to the manufacturer's instructions.

[CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

[WARRANTY]

KURARAY MEDICAL INC. will replace any product that is proved to be defective. KURARAY MEDICAL INC. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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CLEARFIL and PANAVIA are trademarks of KURARAY CO., LTD.
CLEARFIL y PANAVIA son marcas registradas de KURARAY CO., LTD.