CLEARFIL™ AP-X PLT**Intro Kit****ENGLISH INSTRUCTIONS FOR USE****Description**

CLEARFIL AP-X PLT is a light-cured, radiopaque restorative composite resin which provides accurate color-matching, high polishability and excellent physical properties, making it ideal for both anterior and posterior restorations. It is formulated with optimal viscosity assuring easy handling and placement. CLEARFIL AP-X PLT, with the special dispensing system, can be quickly placed directly into the cavity.

CLEARFIL SE BOND is a light-cured bonding system that consists of a self-etching primer and a bonding agent. The primer provides simultaneous treatment of both dentin and enamel.

I. Indications

CLEARFIL AP-X PLT is indicated for the following restorative applications:

1. Class I, II, V restorations of posterior teeth
2. Class III, IV, V restorations of anterior teeth
3. Cervical cavities or root surfaces defects

II. Contraindications

Patients with allergies to methacrylate monomers.

III. Possible Side Effects

The mucous membrane may turn slightly whitish when from contact with the PRIMER due to coagulation.

IV. Incompatibility

Eugenol containing materials for pulp protection and temporary sealing should not be used since they retard the curing process.

V. Precautions**1. Safety precautions**

- Avoid use of the product for patients with a history of hypersensitivity to methacrylate monomers.
- If any hypersensitivity, such as a rash or dermatitis, occurs from with the product, discontinue the use of the product and consult a physician.
- Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers, acetone or acid.
- Use caution to prevent the product from coming in contact with the skin or getting into the eyes. Before using the product, cover the patient's eyes with a towel to protect the patient eyes from possible splashing material.
- If the product comes in contact with human body tissues, take the following actions:
<If the product gets in the eyes>
Immediately wash the eye with copious amounts of water and consult a physician.
<If the product comes in contact with the skin>
Immediately wipe the material off using a cotton pledget (or gauze) moistened with alcohol and then wash with copious amounts of water.
- Use caution to prevent the patient from accidentally swallowing the product.
- Avoid looking directly at the curing light when curing the product.
- Keep the product out of reach of children.

2. Handling and manipulation precautions

- Do not use CLEARFIL AP-X PLT in conjunction with the other composite resin. Mixing materials may cause a change in physical properties, possibly a decrease, from the properties expected.
- When light curing CLEARFIL AP-X PLT, note the light curing depth in this Instructions For Use.
- CLEARFIL AP-X PLT and CLEARFIL SE BOND polymerize by visible light (especially ultraviolet rays). Avoid operating light or natural light (sunlight from windows) and use within 3 minutes of dispensing.
- The emitting tip of the visible light-curing unit should be held as close to and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
- Low light intensity causes poor adhesion and curing. Check the lamp for service life and the visible light tip for contamination. It is advisable to check the curing light intensity at periodic intervals to assure optimal intensity and curing.

3. Storage precautions

- Do not use after the expiration date. Note the expiration date on the outside of the package.
- CLEARFIL AP-X PLT should be stored below 25°C when not being used. CLEARFIL SE BOND should also be refrigerated when not in use, and should be brought to room temperature before using.
- Keep away from extreme heat or direct sunlight.
- Replace the bottle cap as soon as possible after the resin has been dispensed.

VI. System Shades and Components**1. Shades**

CLEARFIL AP-X PLT is available in 8 shades; 6 shades correspond on a best-match basis with VITA Shade Guide plus 2 shades for special requirements. Select an appropriate shade by matching the shade guide to the restoration site. The XL and HO shades offer a strong, lighter white color. The HO shade is translucent white and is suitable for sites requiring extra white color. Standard shades: A2, A3, A3.5, B2, B3
Cervical shades: A4
Translucent shades: XL, HO

2. Kit Components**CLEARFIL AP-X PLT**

10 tips (0.2g/tip) : A2, A3, B2, XL

Note: All 8 shades are available separately.

CLEARFIL SE BOND

PRIMER:.....1.5ml

BOND:.....1.25ml

VII. Clinical Procedures

Use of CLEARFIL SE BOND is especially recommended, although all CLEARFIL bonding agents can be used. Follow the Instructions for use for the specific bonding agent chosen.

1. Cleaning tooth structure

Adequate cavity cleaning is necessary to assure optimal adhesive performance. Make sure the cavity is carefully cleaned.

2. Moisture control

To assure optimal results, avoid blood or saliva contamination of the treatment area. A rubber dam is recommended to keep the tooth dry and clean.

3. Cavity preparations

Remove any infected dentin using CRIES DETECTOR and prepare the cavity in the usual manner.

4. Pulp protection

Any actual or near pulp exposure should be covered with a hard setting calcium hydroxide material. However, usually a cement lining or basing is not necessary. Do not use eugenol materials for pulp protection.

5. Shade selection

Teeth should be cleaned to remove plaque or stains. Before isolating the teeth, select an appropriate shade using the VITA Lumin VACUUM shade guide.

6. Acid etching uncut enamel

If the resin restorative material will extend to uncut enamel, apply K-ETCHANT GEL to the enamel, let it stay for 10 seconds, wash with water, and then dry.

[CAUTION]

The use of PRIMER alone does not adequately condition uncut enamel. Overfilling of resin onto unetched, uncut enamel could cause marginal discoloration.

7. Tooth surface treatment

1) Dispense the necessary amount of PRIMER into the mixing dish well immediately prior to application.

[CAUTION]

PRIMER polymerizes by visible light (especially ultraviolet rays). PRIMER will set to a gel if left under natural light or operating light.

Use the light blocking plate to avoid operating light or natural light (sunlight from windows) and use within 3 minutes of dispensing.

2) Apply PRIMER to the entire cavity wall with a sponge or a disposable brush tip. Leave it in place for 20 seconds. Use caution to prevent saliva or exudate from contacting the treated surfaces for at least 20 seconds.

3) After conditioning the tooth surface for 20 seconds, evaporate the volatile ingredients with a mild oil-free air stream.

[CAUTION]

Avoid pooling of PRIMER. Also do not wash the treated area.

Observe the drying instructions and treatment time to ensure optimum adhesion.

Avoid touching the treated surfaces. If a treated surface is contaminated, wash it with water, dry, or clean with alcohol, and again treat with PRIMER.

8. Bonding

1) Dispense the necessary amount of BOND into a mixing dish well.

[CAUTION]

BOND polymerizes by visible light (especially ultraviolet rays). BOND will set to a gel if left under the operating light or natural light.

Use the light blocking plate to avoid operating light or natural light (sunlight from windows) and use within 3 minutes of dispensing.

2) Apply BOND to the entire surface of the cavity with a sponge or a disposable brush tip.

3) After application, make the bond film as uniform as possible using a gentle oil-free air stream.

[CAUTION]

A strong air blast will scatter the bonding agent, resulting in irregular and poor adhesion.

4) Light-cure the BOND for 10 seconds with a visible light curing unit.

9. Placement and light curing of CLEARFIL AP-X PLT**1) Dispensing**

Place a PLT (Pre-loaded tip) into the dispenser barrel according to the instructions for use for dispenser. (Centrix syringe such as CR EZ or MARK IIIP is recommended.)

The tip may be rotated to gain the proper angle for entrance into the cavity.

Squeeze the dispenser with a slow, steady pressure. Excessive force is not necessary. Remove and discard the tip after use.

2) Placement

Incremental placement and light-curing each increment is strongly recommended especially in deep cavities and in Class II cavities.

3) Curing

Cure the resin with a dental visible light curing unit.

Hold the light tip as close to the resin as possible.

curing depth (mm)

Shades	20sec.	40sec.
A2, A3, B2, XL, HO	2.0	2.5
A3.5, A4, B3	1.5	2.0

[NOTE]

Use of light curing units conforming to ISO/TS 10650:1999 is recommended. If necessary, confirm the curing depth of your unit using method 7.10 of ISO 4049:2000.

10. Finishing

Contour the restoration and adjust the occlusion using a fine diamond point. Polish with silicon rubber points or polishing discs.

[CAUTION]

Federal (U.S.A) law restricts this device to sale by or on the order of a licensed dentist.

[WARRANTY]

KURARAY MEDICAL INC. will replace any product that is shown to be defective.

KURARAY MEDICAL INC. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products.

Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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CLEARFIL is a trade mark of KURARAY CO., LTD.
VITA is a trade mark of VITA Zahnfabrik, H.Rauter GmbH & CO KG Bad Sackingen, Germany.
Las leyes federales (U.S.A) limitan la venta de este producto a dentistas con licencia personas autorizadas por ellos.
CLEARFIL es una marca registrada de KURARAY CO., LTD.
VITA es una marca registrada de VITA Zahnfabrik, H.Rauter GmbH & CO KG Bad Sackingen, Alemania.